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Temporary root canal cement

Composition

Powder Zinc Oxide, Zinc Acetate, Polymethylmethacrylate

Liquid Eugenol

Properties

The therapeutic value of Zinc Oxide-Eugenol cement Pro-Pul-Pan is based on its compatibility with both hard and soft tissues of the mouth. Zinc Oxide-Eugenol cements, like Pro-Pul-Pan, are particularly used upon pulpal tissues, especially when restoring teeth with deep carious lesions. Additionally, it provides a good marginal seal and thermal insulation.

Indications

- > Temporary fillings, thermal insulation, pulp capping
- > Periodontal dressings
- > Root canal sealer
- > Temporary sealing of crowns and inlays

Directions for use

- > Place powder and liquid on a mixing pad at a recommended ratio of 6 to 1 per weight
- > Three spatulation cycles are required to obtain a smooth, creamy mass of the desired consistency
 - Half the powder should be mixed in the first cycle and the remainder in the two subsequent cycles
 - Speed is not essential as setting will only begin 3 or 4 minutes after initiating the mix
 - Mixing on a cool, dry slab induces a slower set. A thin mix also means a slow set
 - Mixing on a warm, damp slab induces a faster set. A thick mix also sets faster
- > The presence of Eugenol will interfere with resin-based fillings

Storage

Store at room temperature between 5°C and 30°C, protected from direct light

Packaging

	REF	REF
	White	Ivory
Portion of 30 g powder + 15 ml liquid	10260	10263
Bottle of 30 g powder	10261	10264
Bottle of 15 ml liquid	10262	10262

> FOR DENTAL PROFESSIONAL USE > MADE IN SWITZERLAND

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